

# The Versatile One !

## Wound Vacuum System



### For the Promotion of Wound Healing

Designed to assist with a variety of wound healing needs. The Versatile 1 Wound Vacuum System® applies a localized negative pressure on the wound bed. This occurs through a high quality piston operation that is assisted by very innovative suction regulator.

- Whisper Quiet and User Friendly.
- Approved indications for use include “may promote wound healing”
- No Annoying Alarms
- Low Levels of Suction (60-80 mmHg) increases patient compliance
- Can be used for multiple applications. Ideal for Military or Triage Use

The Versatile 1™ is the pump for your hospital, institution, or home healthcare agency. It is an economical, robust, and reliable Wound pump that can assist in the promotion of wound healing or wound drainage in many situations. But more than that, Versatile 1™ can be used for virtually any suction need. From intermittent gastric suction to endoscopic use the Versatile 1 is the pump for you.

Versatile 1 Wound Vacuum System is the only approach that offers so many options to clinicians. The product can be used for the promotion of wound healing. It can also be used for many complex wound drainage applications. This is the ideal pump for Advanced Practitioners.

### Wound Vacuum System

6965 El Camino Real Suite 105-602  
La Costa, CA 92009  
United States of America  
Info@blueskymedical.com

760-727-1477  
760-727-0236 (FAX)



# Versatile 1 Wound Vacuum System™

## *Its Simply Better!*

### *Chariker-Jeter® Wound Sealing Kit*

In 1985 Katherine Jeter, Ed.D, ET came up with a unique combination of products to apply negative pressure to the wound bed. Dr. Mark Chariker helped to put a clinical study together and further this art. They knew that application of negative pressure to the wound bed had extraordinary effects. They and others published this information for the general good of patients and to advance the art of wound care in 1989.

The Chariker-Jeter Kit Wound Sealing Kit is an exact embodiment of this kit. Classical components are used. The kit has a comprehensive set of instructions. There is extensive documentation plus color step by step instructions and protocols on our website. BlueSky believes in education for healthcare professionals with its products.

The kit is available in several forms:

500.7777 Chariker-Jeter Wound Sealing Kit

500.2005 Chariker-Jeter Basic Wound Sealing Kit

500.3010 Chariker-Jeter Basic Plus Kit

500. 2045 Chariker-Jeter Jade Kit

#### **Why is It Simply Better?**

- No reports of tissue tearing
- Kits can stay on longer
- Simple approach (Build up)
- Cost effective
- Lower levels of pressure
- No irritating alarms
- Additional applications
- Anti-microbial fill
- Treats wounds large and small
- More comfortable drape
- Periodic use available
- Ideal for military or triage use
- Tunneling Wound Management
- Fistula Management and Drainage


#### **PH Strips**



BlueSky PH Strips accurately chart Wound Fluid PH. Wound treated with


negative pressure should become more alkaline. Item #700.0806

#### **Suction Canister**



The 800ml disposable suction canister is the ideal accessory. Item# 700.9005

#### **Bacteria/Overflow Filter**



Provides protection to the pump against overflow. Also is a bacteria filter. Item #500.3070

#### **Wound Vacuum System**

6965 El Camino Real Suite 105-602  
La Costa, CA 92009  
United States of America  
Info@blueskymedical.com

760-727-1477  
760-727-0236 (FAX)

